

产品名称: ApoA4 (5G1) Rabbit Monoclonal Antibody
产品货号: AMRe07019



产品概述 (Summary)

产品名称 (Production Name)	ApoA4 (5G1) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOA4
别名 (Alternative Names)	Apo AIV; APOA 4; ApoA IV; APOA4; Apolipoprotein A IV; Apolipoprotein A4;
基因 ID (Gene ID)	337.0
蛋白 ID (SwissProt ID)	P06727.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)

产品名称: ApoA4 (5G1) Rabbit Monoclonal Antibody
产品货号: AMRe07019




APOA4 (apolipoprotein A-IV) is a component of HDL and chylomicrons. Its primary site of synthesis is the intestine, in association with lymph chylomicron particles. May have a role in chylomicrons and VLDL secretion and catabolism. Required for efficient activation of lipoprotein lipase by ApoC-II; potent activator of LCAT. ApoA-IV is a major component of HDL and chylomicrons.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of ApoA4 expression in Human plasma lysate.

注意事项 (Note)

For research use only .