

产品名称: ApoD (18D19) Rabbit Monoclonal Antibody
产品货号: AMRe07033



产品概述 (Summary)

产品名称 (Production Name)	ApoD (18D19) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOD
别名 (Alternative Names)	ApoD; ApolipoproteinD;
基因 ID (Gene ID)	347.0
蛋白 ID (SwissProt ID)	P05090.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	21kDa

研究背景 (Background)

产品名称: ApoD (18D19) Rabbit Monoclonal Antibody
产品货号: AMRe07033




APOD occurs in the macromolecular complex with lecithin-cholesterol acyltransferase. It is probably involved in the transport and binding of bilin. Appears to be able to transport a variety of ligands in a number of different contexts. APOD occurs in the macromolecular complex with lecithin- cholesterol acyltransferase. It is probably involved in the transport and binding of bilin. Appears to be able to transport a variety of ligands in a number of different contexts.

研究领域 (Research Area)

Neuroscience; Neurology process; Neurodegenerative disease; Alzheimer's disease; Cardiovascular; Lipids / Lipoproteins; Lipoproteins/Apolipoproteins; Signal Transduction; Metabolism Lipid metabolism; Cancer; Tumor biomarkers; Neurodegenerative disease

图片 (Image Data)



Western blot analysis of ApoD expression in Human plasma lysate.

注意事项 (Note)

For research use only .