

产品名称: ARF6 (19O13) Rabbit Monoclonal Antibody
产品货号: AMRe07102



产品概述 (Summary)

产品名称 (Production Name)	ARF6 (19O13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARF6
别名 (Alternative Names)	ADP ribosylation factor 6; ARF6; Small GTP binding protein;
基因 ID (Gene ID)	382.0
蛋白 ID (SwissProt ID)	P62330.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	20kDa

研究背景 (Background)

产品名称: ARF6 (19O13) Rabbit Monoclonal Antibody
产品货号: AMRe07102




ADP-ribosylation factor (Arf) proteins are low molecular weight GTP binding proteins that belong to the Ras GTPase superfamily. May modulate vesicle budding and uncoating within the Golgi apparatus. Functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in the regulation of dendritic spine development (By similarity). GTP-binding protein involved in protein trafficking that regulates endocytic recycling and cytoskeleton remodeling (PubMed:11266366, PubMed:21170023, PubMed:16737952, PubMed:7589240, PubMed:18400762). Required for normal completion of mitotic cytokinesis (By similarity). Plays a role in the reorganization of the actin cytoskeleton and the formation of stress fibers (By similarity). Involved in the regulation of dendritic spine development, contributing to the regulation of dendritic branching and filopodia extension (PubMed:14978216). Plays an important role in membrane trafficking, during junctional remodeling and epithelial polarization. Regulates surface levels of adherens junction proteins such as CDH1 (By similarity). Required for NTRK1 sorting to the recycling pathway from early endosomes (By similarity).

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of ARF6 expression in 293T cell lysate.

注意事项 (Note)

For research use only .