

产品名称: BAG2 (17W11) Rabbit Monoclonal Antibody
产品货号: AMRe07435



产品概述 (Summary)

产品名称 (Production Name)	BAG2 (17W11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAG2
别名 (Alternative Names)	BAG 2; Bag2;
基因 ID (Gene ID)	9532.0
蛋白 ID (SwissProt ID)	O95816.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)

产品名称: BAG2 (17W11) Rabbit Monoclonal Antibody
产品货号: AMRe07435




Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Co-chaperone for HSP70 and HSC70 chaperone proteins. Acts as a nucleotide-exchange factor (NEF) promoting the release of ADP from the HSP70 and HSC70 proteins thereby triggering client/substrate protein release (PubMed:24318877, PubMed:9873016).

研究领域 (Research Area)

Cell Biology; Apoptosis; Intracellular; Bcl2 Family; Signal Transduction; Protein Trafficking; Chaperones; Cancer; Cell Death

图片 (Image Data)



Western blot analysis of BAG2 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .