

产品名称: BAG2 (9Z17) Rabbit Monoclonal Antibody
产品货号: AMRe07436



产品概述 (Summary)

产品名称 (Production Name)	BAG2 (9Z17) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAG2
别名 (Alternative Names)	Bag2;
基因 ID (Gene ID)	9532.0
蛋白 ID (SwissProt ID)	O95816.A synthetic peptide of human BAG2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:20-1:100,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:20-1:50,IF-P 1:20-1:50
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)

产品名称: BAG2 (9Z17) Rabbit Monoclonal Antibody
产品货号: AMRe07436




Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Co-chaperone for HSP70 and HSC70 chaperone proteins. Acts as a nucleotide-exchange factor (NEF) promoting the release of ADP from the HSP70 and HSC70 proteins thereby triggering client/substrate protein release (PubMed:24318877, PubMed:9873016).

研究领域 (Research Area)

Cell Biology; Apoptosis; Intracellular; Bcl2 Family; Signal Transduction; Protein Trafficking; Chaperones; Cancer; Cell Death

图片 (Image Data)



Western blot detection of BAG2 in Jurkat, HeLa cell lysates using BAG2 antibody (1:1000 diluted).

注意事项 (Note)

For research use only .