

**产品名称:** Bak (16T6) Rabbit Monoclonal Antibody

**产品货号:** AMRe07456

### 产品概述 (Summary)

|                   |                                                                                                                             |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 描述(Description)   | Recombinant rabbit monoclonal antibody                                                                                      |
| 宿主(Host)          | Rabbit                                                                                                                      |
| 应用(Application)   | WB,IHC,ICC/IF,FC,IP                                                                                                         |
| 种属反应(Reactivity)  | Human                                                                                                                       |
| 偶联物(Conjugation)  | Unconjugated                                                                                                                |
| 修饰(Modification)  | Unmodified                                                                                                                  |
| 同种型(Isotype)      | IgG                                                                                                                         |
| 克隆(Clonality)     | Monoclonal                                                                                                                  |
| 剂型(Form)          | Liquid                                                                                                                      |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |
| 储存溶液(Buffer)      | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| 纯化 (Purification) | Affinity purification                                                                                                       |

### 产品应用(Application)

|                       |                                                                              |
|-----------------------|------------------------------------------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:2000-1:10000,IHC 1:100-1:200,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:20-1:50 |
| 分子量(Molecular Weight) | 23kDa                                                                        |

### 抗原信息(Antigen Information)

|                       |                                  |
|-----------------------|----------------------------------|
| 基因名(Gene Name)        | BAK1                             |
| 别名(Alternative Names) | BAK1;BAK;BAK-LIKE;BCL2L7;CDN1;   |
| 基因 ID(Gene ID)        | 578.0                            |
| SwissProt ID          | Q16611                           |
| 免疫原(Immunogen)        | A synthetic peptide of human Bak |

### 研究背景 (Background)

Bak is a proapoptotic member of the Bcl-2 family. This protein is located on the outer membrane of mitochondria and is an essential component for transduction of apoptotic signals through the mitochondrial pathway. Upon apoptotic stimulation, an upstream stimulator like truncated BID (tBID) induces conformational changes in Bak to form oligomer channels in the

mitochondrial membrane for cytochrome c release. The release of cytochrome c to the cytosol activates the caspase-9 pathway and eventually leads to cell death. Plays a role in the mitochondrial apoptotic process. Upon arrival of cell death signals, promotes mitochondrial outer membrane (MOM) permeabilization by oligomerizing to form pores within the MOM. This releases apoptogenic factors into the cytosol, including cytochrome c, promoting the activation of caspase 9 which in turn processes and activates the effector caspases.


### 研究领域 (Research Area)

Cell Biology; Apoptosis; Intracellular; Bcl2 Family; Oligo verified

### 注意事项 (Note)

For research use only.

### 图片 (Image Data)



Western blot detection of Bak in K562, HeLa cell lysates using Bak antibody (1:1000 diluted).

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>