

产品名称: Bid (6Y3) Rabbit Monoclonal Antibody
产品货号: AMRe07557



产品概述 (Summary)

产品名称 (Production Name)	Bid (6Y3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BID
别名 (Alternative Names)	BID isoform ES(1b); BID isoform L(2); BID isoform Si6; FP497; p11 BID; p13 BID; p15 BID; p22 BID;
基因 ID (Gene ID)	637.0
蛋白 ID (SwissProt ID)	P55957.A synthetic peptide of human Bid

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:20-1:100,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:20-1:50,IF-P 1:20-1:50
蛋白分子量 (Molecular Weight)	22kDa

产品名称: Bid (6Y3) Rabbit Monoclonal Antibody
产品货号: AMRe07557




研究背景 (Background)

The major proteolytic product p15 BID allows the release of cytochrome c (By similarity). Isoform 1, isoform 2 and isoform 4 induce ICE-like proteases and apoptosis. Isoform 3 does not induce apoptosis. Counters the protective effect of Bcl-2. Induces caspases and apoptosis (PubMed: <http://www.uniprot.org/citations/14583606> target="_blank">14583606). Counters the protective effect of BCL2 (By similarity).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of Bid in HeLa,A549,HL-60 using Bid antibody(1:1000 diluted).

注意事项 (Note)

For research use only .