

产品名称: C3orf38 (11N12) Rabbit Monoclonal Antibody
产品货号: AMRe07748



产品概述 (Summary)

产品名称 (Production Name)	C3orf38 (11N12) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C3orf38
别名 (Alternative Names)	C3orf38;
基因 ID (Gene ID)	285237.0
蛋白 ID (SwissProt ID)	Q5JPI3.Recombinant protein of human C3orf38

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	38kDa

研究背景 (Background)

May be involved in apoptosis regulation. May be involved in apoptosis regulation.


产品名称: C3orf38 (11N12) Rabbit Monoclonal Antibody
产品货号: AMRe07748



研究领域 (Research Area)

Cell Biology; Apoptosis; Intracellular; Associated Proteins

图片 (Image Data)



Western blot detection of C3orf38 in K562, HeLa cell lysates using C3orf38 antibody (1:1000 diluted).

注意事项 (Note)

For research use only .