

产品名称: Cadherin 16 (15J2) Rabbit Monoclonal Antibody
产品货号: AMRe07826



产品概述 (Summary)

产品名称 (Production Name)	Cadherin 16 (15J2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDH16
别名 (Alternative Names)	Cadherin 16; Cadherin kidney; CDH16; Kidney cadherin; Ksp cadherin;
基因 ID (Gene ID)	1014.0
蛋白 ID (SwissProt ID)	O75309.Recombinant protein of human Cadherin 16

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:5000,IP 1:50-1:100,IF-P 1:500-1:5000
蛋白分子量 (Molecular Weight)	90kDa

研究背景 (Background)

Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic

产品名称: Cadherin 16 (15J2) Rabbit Monoclonal Antibody
产品货号: AMRe07826




manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

研究领域 (Research Area)

Signal Transduction; Cytoskeleton / ECM; Cell Adhesion; Cadherins

图片 (Image Data)



Western blot analysis of extracts from Mouse kidney tissue using RM6027 at 1:1000.

注意事项 (Note)

For research use only .