

产品名称: CAP1 (14W1) Rabbit Monoclonal Antibody
产品货号: AMRe07909



产品概述 (Summary)

产品名称 (Production Name)	CAP1 (14W1) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAP1
别名 (Alternative Names)	CAP1; CAP1 PEN;
基因 ID (Gene ID)	10487.0
蛋白 ID (SwissProt ID)	Q01518.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200,ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	52kDa

研究背景 (Background)

产品名称: CAP1 (14W1) Rabbit Monoclonal Antibody
产品货号: AMRe07909




Directly regulates filament dynamics and has been implicated in a number of complex developmental and morphological processes, including mRNA localization and the establishment of cell polarity. Directly regulates filament dynamics and has been implicated in a number of complex developmental and morphological processes, including mRNA localization and the establishment of cell polarity.

研究领域 (Research Area)

Signal Transduction; Cytoskeleton / ECM; Extracellular Matrix; Structures; Focal Adhesions; Epigenetics and Nuclear Signaling; DNA / RNA; RNA Processing

图片 (Image Data)



Western blot analysis of CAP1 expression in Raji cell lysate.

注意事项 (Note)

For research use only .