

产品名称: CCL23 (3Y14) Rabbit Monoclonal Antibody
产品货号: AMRe08139



产品概述 (Summary)

产品名称 (Production Name)	CCL23 (3Y14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCL23
别名 (Alternative Names)	CCL23; CKB8; Hmrp-2a; MIP3; MPIF1; SCYA23;
基因 ID (Gene ID)	6368.0
蛋白 ID (SwissProt ID)	P55773.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,IF-P 1:50-1:100
蛋白分子量 (Molecular Weight)	13kDa

研究背景 (Background)


产品名称: CCL23 (3Y14) Rabbit Monoclonal Antibody
产品货号: AMRe08139



Shows chemotactic activity for monocytes, resting T-lymphocytes, and neutrophils, but not for activated lymphocytes. Inhibits proliferation of myeloid progenitor cells in colony formation assays. This protein can bind heparin. Binds CCR1. CCL23(19-99), CCL23(22-99), CCL23(27-99), CCL23(30-99) are more potent chemoattractants than the small-inducible cytokine A23. Shows chemotactic activity for monocytes, resting T- lymphocytes, and neutrophils, but not for activated lymphocytes. Inhibits proliferation of myeloid progenitor cells in colony formation assays. This protein can bind heparin. Binds CCR1. CCL23(19-99), CCL23(22-99), CCL23(27-99), CCL23(30-99) are more potent chemoattractants than the small-inducible cytokine A23.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CCL23 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .