

产品名称: CCR8 (11B19) Rabbit Monoclonal Antibody
产品货号: AMRe08161



产品概述 (Summary)

产品名称 (Production Name)	CCR8 (11B19) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCR8
别名 (Alternative Names)	CY6; TER1; CCR-8; CKRL1; CDw198; CMKBR8; GPRCY6; CMKBRL2; CC-CCR-8; TER1;
基因 ID (Gene ID)	1237.0
蛋白 ID (SwissProt ID)	P51685.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200, FC 1:500-1:2000
蛋白分子量 (Molecular Weight)	41kDa

产品名称: CCR8 (11B19) Rabbit Monoclonal Antibody
产品货号: AMRe08161




研究背景 (Background)

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection. Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection.

研究领域 (Research Area)

Immunology; Innate Immunity; Chemokines; Beta Chemokine Rec. (CCR); Immune System Diseases; Antiviral Signaling; HIV-related; Microbiology; Interspecies Interaction; Host Virus Interaction

图片 (Image Data)



Western blot analysis of CCR8 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .