

产品名称: CD19 (16P1) Rabbit Monoclonal Antibody  
产品货号: AMRe08254



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | CD19 (16P1) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,FC                       |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                                                     |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 纯化方式 (Purification) | Affinity purification                                                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                             |
|------------------------|-------------------------------------------------------------|
| 基因名 (Gene Name)        | CD19                                                        |
| 别名 (Alternative Names) | B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; |
| 基因 ID (Gene ID)        | 930.0                                                       |
| 蛋白 ID (SwissProt ID)   | P15391.                                                     |

## 产品应用 (Application)

|                          |                                                                      |
|--------------------------|----------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:200-1:1000,FC 1:500-1:1000 |
| 蛋白分子量 (Molecular Weight) | 61kDa                                                                |

## 研究背景 (Background)

产品名称: CD19 (16P1) Rabbit Monoclonal Antibody  
产品货号: AMRe08254




CD19 is a transmembrane glycoprotein that contains two extracellular immunoglobulin-like domains. Assembles with the antigen receptor of B lymphocytes in order to decrease the threshold for antigen receptor-dependent stimulation. Functions as coreceptor for the B-cell antigen receptor complex (BCR) on B-lymphocytes. Decreases the threshold for activation of downstream signaling pathways and for triggering B-cell responses to antigens (PubMed:<a href="http://www.uniprot.org/citations/2463100" target="\_blank">2463100</a>, PubMed:<a href="http://www.uniprot.org/citations/1373518" target="\_blank">1373518</a>, PubMed:<a href="http://www.uniprot.org/citations/16672701" target="\_blank">16672701</a>). Activates signaling pathways that lead to the activation of phosphatidylinositol 3-kinase and the mobilization of intracellular Ca(2+) stores (PubMed:<a href="http://www.uniprot.org/citations/9382888" target="\_blank">9382888</a>, PubMed:<a href="http://www.uniprot.org/citations/9317126" target="\_blank">9317126</a>, PubMed:<a href="http://www.uniprot.org/citations/12387743" target="\_blank">12387743</a>, PubMed:<a href="http://www.uniprot.org/citations/16672701" target="\_blank">16672701</a>). Is not required for early steps during B cell differentiation in the blood marrow (PubMed:<a href="http://www.uniprot.org/citations/9317126" target="\_blank">9317126</a>). Required for normal differentiation of B-1 cells (By similarity). Required for normal B cell differentiation and proliferation in response to antigen challenges (PubMed:<a href="http://www.uniprot.org/citations/2463100" target="\_blank">2463100</a>, PubMed:<a href="http://www.uniprot.org/citations/1373518" target="\_blank">1373518</a>). Required for normal levels of serum immunoglobulins, and for production of high-affinity antibodies in response to antigen challenge (PubMed:<a href="http://www.uniprot.org/citations/9317126" target="\_blank">9317126</a>, PubMed:<a href="http://www.uniprot.org/citations/12387743" target="\_blank">12387743</a>, PubMed:<a href="http://www.uniprot.org/citations/16672701" target="\_blank">16672701</a>).

## 研究领域 (Research Area)

Immunology

## 图片 (Image Data)



Western blot analysis of CD19 expression in Ramos cell lysate.

产品名称: CD19 (16P1) Rabbit Monoclonal Antibody  
产品货号: AMRe08254

---



**注意事项 (Note)**

For research use only .