

产品名称: CD1a (11L3) Rabbit Monoclonal Antibody
产品货号: AMRe08258



产品概述 (Summary)

产品名称 (Production Name)	CD1a (11L3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD1A
别名 (Alternative Names)	CD 1a; CD1; CD1a; epidermal dendritic cell marker CD1a; FCB6; HTA1; hTa1 thymocyte antigen; R4; T-cell surface antigen T6/Leu-6; T6;
基因 ID (Gene ID)	909.0
蛋白 ID (SwissProt ID)	P06126.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,IF-P 1:200-1:1000
蛋白分子量 (Molecular Weight)	37kDa

产品名称: CD1a (11L3) Rabbit Monoclonal Antibody
产品货号: AMRe08258




研究背景 (Background)

Antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells. Antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD1a expression in human fetal thymus lysate.

注意事项 (Note)

For research use only .