

产品名称: CD23 (17L19) Rabbit Monoclonal Antibody
产品货号: AMRe08292



产品概述 (Summary)

产品名称 (Production Name)	CD23 (17L19) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FCER2
别名 (Alternative Names)	CD23; FCE2; CD23A; IGEBF; CLEC4J; FCER2;
基因 ID (Gene ID)	2208.0
蛋白 ID (SwissProt ID)	P06734.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,FC 1:20-1:50,IF-P 1:100-1:200
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)

产品名称: CD23 (17L19) Rabbit Monoclonal Antibody
产品货号: AMRe08292




This receptor has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen). Low-affinity receptor for immunoglobulin E (IgE) and CR2/CD21. Has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD23 expression in Raji cell lysate.

注意事项 (Note)

For research use only .