

产品名称: CD32 (12Q18) Rabbit Monoclonal Antibody
产品货号: AMRe08359



产品概述 (Summary)

产品名称 (Production Name)	CD32 (12Q18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FCGR2A
别名 (Alternative Names)	CD32; CD32A; CD32B; CD32C; FCG2; FcGR; FCGR2; FCGR2A; FCGR2B; FCGR2C; Fcgr3a; FCRII; IGFR2;
基因 ID (Gene ID)	2212.0
蛋白 ID (SwissProt ID)	P12318.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,IF-P 1:50-1:100
蛋白分子量 (Molecular Weight)	35kDa

产品名称: CD32 (12Q18) Rabbit Monoclonal Antibody
产品货号: AMRe08359




研究背景 (Background)

Binds to the Fc region of immunoglobulins gamma. Low affinity receptor. By binding to IgG it initiates cellular responses against pathogens and soluble antigens. Binds to the Fc region of immunoglobulins gamma. Low affinity receptor. By binding to IgG it initiates cellular responses against pathogens and soluble antigens. Promotes phagocytosis of opsonized antigens.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of CD32 expression in K562 cell lysate.

注意事项 (Note)

For research use only .