

产品名称: CD97 (16U13) Rabbit Monoclonal Antibody
产品货号: AMRe08491



产品概述 (Summary)

产品名称 (Production Name)	CD97 (16U13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADGRE5
别名 (Alternative Names)	ADGRE5; Adhesion G protein-coupled receptor E5; CD97 antigen; Leukocyte antigen CD97; Seven span transmembrane protein; TM7LN1;
基因 ID (Gene ID)	976.0
蛋白 ID (SwissProt ID)	P48960.A synthetic peptide of human CD97

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	92kDa

研究背景 (Background)


产品名称: CD97 (16U13) Rabbit Monoclonal Antibody
产品货号: AMRe08491



Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration. Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Jurkat cells using RM5542 at 1:1000.

注意事项 (Note)

For research use only .