

产品名称: CDKN2C (10M18) Rabbit Monoclonal Antibody
产品货号: AMRe08579



产品概述 (Summary)

产品名称 (Production Name)	CDKN2C (10M18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDKN2C
别名 (Alternative Names)	Cdkn2c; CDKN6; INK4C; p18; p18 inhibits CDK4; p18 INK4c; p18 INK6; p18-INK4c; p18-INK6;
基因 ID (Gene ID)	1031.0
蛋白 ID (SwissProt ID)	P42773.Recombinant protein of human p18 INK4c

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	18kDa

研究背景 (Background)

产品名称: CDKN2C (10M18) Rabbit Monoclonal Antibody
产品货号: AMRe08579




Interacts strongly with CDK6, weakly with CDK4. Inhibits cell growth and proliferation with a correlated dependence on endogenous retinoblastoma protein RB. Interacts strongly with CDK6, weakly with CDK4. Inhibits cell growth and proliferation with a correlated dependence on endogenous retinoblastoma protein RB.

研究领域 (Research Area)

Cell_Cycle_G1S;Cell_Cycle_G2M_DNA;

图片 (Image Data)



Western blot detection of p18 INK4c/CDKN2C in C6 cell lysates using p18 INK4c/CDKN2C antibody(1:500 diluted).

注意事项 (Note)

For research use only .