

产品名称: CECR5 (17M19) Rabbit Monoclonal Antibody
产品货号: AMRe08624



产品概述 (Summary)

产品名称 (Production Name)	CECR5 (17M19) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HDHD5
别名 (Alternative Names)	Cecr5;
基因 ID (Gene ID)	27440.0
蛋白 ID (SwissProt ID)	Q9BXW7.A synthetic peptide of human CECR5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	46kDa

研究背景 (Background)

Candidate gene for the Cat Eye Syndrome (CES), a developmental disorder associated with the duplication of a 2 Mb region

产品名称: CECR5 (17M19) Rabbit Monoclonal Antibody
产品货号: AMRe08624




of 22q11.2. Duplication usually takes in the form of a supernumerary bisatellited isodicentric chromosome, resulting in four copies of the region (represents an inv dup(22)(q11)). CES is characterized clinically by the combination of coloboma of the iris and anal atresia with fistula, downslanting palpebral fissures, preauricular tags and/or pits, frequent occurrence of heart and renal malformations, and normal or near-normal mental development.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of CECR5 in K562, HeLa cell lysates using CECR5 antibody (1:1000 diluted).

注意事项 (Note)

For research use only.