

产品名称: CEP43 (7K4) Rabbit Monoclonal Antibody
产品货号: AMRe08665



产品概述 (Summary)

产品名称 (Production Name)	CEP43 (7K4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEP43
别名 (Alternative Names)	FGFR1OP; FOP;
基因 ID (Gene ID)	11116.0
蛋白 ID (SwissProt ID)	O95684.A synthetic peptide of human FGFR1 Oncogene Partner

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:20-1:50, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	43kDa

研究背景 (Background)

Required for anchoring microtubules to the centrosomes. Required for anchoring microtubules to the centrosomes


产品名称: CEP43 (7K4) Rabbit Monoclonal Antibody
产品货号: AMRe08665



(PubMed:16314388, PubMed:28659385). Required for ciliation (PubMed:28625565, PubMed:28659385).

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of FGFR1 Oncogene Partner in A549,HL-60,U251,U87-MG cell lysates using FGFR1 Oncogene Partner antibody(1:1000 diluted).

注意事项 (Note)

For research use only .