

产品名称: CHD3 (5M11) Rabbit Monoclonal Antibody
产品货号: AMRe08745



产品概述 (Summary)

产品名称 (Production Name)	CHD3 (5M11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHD3
别名 (Alternative Names)	Chd3; Chd7; hZFH; Mi2 ALPHA; Prp7; Prp9 1; ZFH;
基因 ID (Gene ID)	1107.0
蛋白 ID (SwissProt ID)	Q12873.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:100-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	227kDa

研究背景 (Background)

产品名称: CHD3 (5M11) Rabbit Monoclonal Antibody
产品货号: AMRe08745




Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity. Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones (PubMed:9804427, PubMed:30397230). Involved in transcriptional repression as part of the NuRD complex (PubMed:27068747). Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity (PubMed:17626165).

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CHD3 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .