

产品名称: Clathrin heavy chain (18K9) Rabbit Monoclonal Antibody
产品货号: AMRe08889



产品概述 (Summary)

产品名称 (Production Name)	Clathrin heavy chain (18K9) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	CHC17; Clathrin heavy chain; Clathrin heavy chain like 1; CLH17; CLH22; CLTC; CLTCL; CLTD;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q00610/P53675.A synthetic peptide of human Clathrin heavy chain

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:100,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	192kDa


产品名称: Clathrin heavy chain (18K9) Rabbit Monoclonal Antibody
产品货号: AMRe08889

研究背景 (Background)

Clathrin heavy chain proteins are composed of several functional domains, including a carboxy-terminal region that permits interaction with other heavy chain proteins within a triskelion, and a globular amino-terminal region that associates with other vesicle proteins. Adaptor proteins, such as AP2, epsin and EPS15, are responsible for the recruitment of vesicle proteins to sites of pit formation and the assembly of the clathrin-coated vesicle.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Clathrin heavy chain in HeLa, A549, HL-60, U251, U87-MG cell lysates using Clathrin heavy chain antibody (1:1000 diluted).

注意事项 (Note)

For research use only .