

产品名称: Claudin 11 (15R12) Rabbit Monoclonal Antibody
产品货号: AMRe08891



产品概述 (Summary)

产品名称 (Production Name)	Claudin 11 (15R12) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLDN11
别名 (Alternative Names)	Claudin 11; CLD11_HUMAN; Cldn11; OSP; OTM;
基因 ID (Gene ID)	5010.0
蛋白 ID (SwissProt ID)	O75508.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:100-1:500
蛋白分子量 (Molecular Weight)	22kDa

研究背景 (Background)


产品名称: Claudin 11 (15R12) Rabbit Monoclonal Antibody
产品货号: AMRe08891



Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Oligodendrocyte Specific Protein expression in SH-SY5Y cell lysate.

注意事项 (Note)

For research use only .