

产品名称: Collagen I (7F14) Rabbit Monoclonal Antibody
产品货号: AMRe09207



产品概述 (Summary)

产品名称 (Production Name)	Collagen I (7F14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	Collagen alpha-1(I) chain, Alpha-1 type I collagen, COL1A1;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P02452/P08123.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	139kDa

研究背景 (Background)


产品名称: Collagen I (7F14) Rabbit Monoclonal Antibody
产品货号: AMRe09207



This gene encodes one of the chains for type I collagen, the fibrillar collagen found in most connective tissues. Mutations in this gene are associated with osteogenesis imperfecta, Ehlers-Danlos syndrome, idiopathic osteoporosis, and atypical Marfan syndrome.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Collagen I in Human stomach tissue lysate.

注意事项 (Note)

For research use only .