

产品名称: Collagen XVII (5G9) Rabbit Monoclonal Antibody
产品货号: AMRe09220



产品概述 (Summary)

产品名称 (Production Name)	Collagen XVII (5G9) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	COL17A1
别名 (Alternative Names)	BP180; BPA 2; BPAG2; COL17A1; Collagen 17; CollagenXVII; LAD1
基因 ID (Gene ID)	1308.0
蛋白 ID (SwissProt ID)	Q9UMD9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:200,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	150kDa

研究背景 (Background)


产品名称: Collagen XVII (5G9) Rabbit Monoclonal Antibody
产品货号: AMRe09220



May play a role in the integrity of hemidesmosome and the attachment of basal keratinocytes to the underlying basement membrane. The 120 kDa linear IgA disease antigen is an anchoring filament component involved in dermal-epidermal cohesion. Is the target of linear IgA bullous dermatosis autoantibodies. May play a role in the integrity of hemidesmosome and the attachment of basal keratinocytes to the underlying basement membrane.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Collagen XVII expression in Human skin lysate.

注意事项 (Note)

For research use only .