

产品名称: CRABP2 (12M4) Rabbit Monoclonal Antibody
产品货号: AMRe09361



产品概述 (Summary)

产品名称 (Production Name)	CRABP2 (12M4) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRABP2
别名 (Alternative Names)	Cellular retinoic acid binding protein 2; Cellular retinoic acid binding protein II ; Crabp2; CRABPII; RBP6;
基因 ID (Gene ID)	1382.0
蛋白 ID (SwissProt ID)	P29373.Recombinant protein of human CRABP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	16kDa

研究背景 (Background)

产品名称: CRABP2 (12M4) Rabbit Monoclonal Antibody
产品货号: AMRe09361




Transports retinoic acid to the nucleus. Regulates the access of retinoic acid to the nuclear retinoic acid receptors.
Transports retinoic acid to the nucleus. Regulates the access of retinoic acid to the nuclear retinoic acid receptors.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of CRABP2 in HeLa cell lysates using CRABP2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .