

产品名称: CSDE1 (16X6) Rabbit Monoclonal Antibody
产品货号: AMRe09447



产品概述 (Summary)

产品名称 (Production Name)	CSDE1 (16X6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CSDE1
别名 (Alternative Names)	Csde1; NRU; UNR;
基因 ID (Gene ID)	7812.0
蛋白 ID (SwissProt ID)	O75534.A synthetic peptide of human CSDE1/NRU

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	89kDa

研究背景 (Background)

RNA-binding protein. Required for internal initiation of translation of human rhinovirus RNA. May be involved in

产品名称: CSDE1 (16X6) Rabbit Monoclonal Antibody
产品货号: AMRe09447




translationally coupled mRNA turnover. RNA-binding protein involved in translationally coupled mRNA turnover (PubMed:11051545, PubMed:15314026). Implicated with other RNA- binding proteins in the cytoplasmic deadenylation/translational and decay interplay of the FOS mRNA mediated by the major coding-region determinant of instability (mCRD) domain (PubMed:11051545, PubMed:15314026). Required for efficient formation of stress granules (PubMed:29395067).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of CSDE1/NRU in A549,HL-60,U2OS using CSDE1/NRU antibody(1:1000 diluted)

注意事项 (Note)

For research use only .