

产品名称: CTLA4 (CD152) (14H2) Rabbit Monoclonal Antibody
产品货号: AMRe09507



产品概述 (Summary)

产品名称 (Production Name)	CTLA4 (CD152) (14H2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CTLA4
别名 (Alternative Names)	CTLA4;CD;CD152;CELIAC3;CTLA-4;GRD4;GSE;ICOS;IDDM12;
基因 ID (Gene ID)	1493.0
蛋白 ID (SwissProt ID)	P16410.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,FC 1:200-1:500,IP 1:20-1:50,IF-P 1:200-1:500
蛋白分子量 (Molecular Weight)	25kDa

产品名称: CTLA4 (CD152) (14H2) Rabbit Monoclonal Antibody
产品货号: AMRe09507




研究背景 (Background)

Inhibitory receptor acting as a major negative regulator of T-cell responses. The affinity of CTLA4 for its natural B7 family ligands, CD80 and CD86, is considerably stronger than the affinity of their cognate stimulatory coreceptor CD28. Inhibitory receptor acting as a major negative regulator of T-cell responses. The affinity of CTLA4 for its natural B7 family ligands, CD80 and CD86, is considerably stronger than the affinity of their cognate stimulatory coreceptor CD28.

研究领域 (Research Area)

Cell adhesion molecules (CAMs);T_Cell_Receptor;Autoimmune thyroid disease;

图片 (Image Data)



Western blot analysis of CTLA4 (CD152) expression in NIH/3T3 cell lysate.

注意事项 (Note)

For research use only .