

产品名称: Cystatin B (7U18) Rabbit Monoclonal Antibody  
产品货号: AMRe09692



## 产品概述 (Summary)

|                        |                                              |
|------------------------|----------------------------------------------|
| 产品名称 (Production Name) | Cystatin B (7U18) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody       |
| 宿主 (Host)              | Rabbit                                       |
| 应用 (Application)       | WB,IHC,IF-P                                  |
| 种属反应性 (Reactivity)     | Human,Rat                                    |

## 产品性能 (Performance)

|                     |                                                                                                                             |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                                                  |
| 形式 (Form)           | Liquid                                                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| 纯化方式 (Purification) | Affinity purification                                                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CSTB                                                                          |
| 别名 (Alternative Names) | CHROW21; CST6; CSTB; Cystatin B (stefin B); CYTB EPM1; EPM1A; PME; STFB; ULD; |
| 基因 ID (Gene ID)        | 1476.0                                                                        |
| 蛋白 ID (SwissProt ID)   | P04080.A synthetic peptide of human Cystatin-B                                |

## 产品应用 (Application)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,IHC 1:100-1:200,IF-P 1:100-1:200 |
| 蛋白分子量 (Molecular Weight) | 11kDa                                             |

## 研究背景 (Background)

**产品名称: Cystatin B (7U18) Rabbit Monoclonal Antibody**  
**产品货号: AMRe09692**




This is an intracellular thiol proteinase inhibitor. Tightly binding reversible inhibitor of cathepsins L, H and B. This is an intracellular thiol proteinase inhibitor. Tightly binding reversible inhibitor of cathepsins L, H and B.

## 研究领域 (Research Area)

Cell Biology

## 图片 (Image Data)



Western blot detection of Cystatin-B in HeLa,A549,HL-60 using Cystatin-B antibody(1:1000 diluted)

## 注意事项 (Note)

For research use only .