

产品名称: Cytochrome P450 Reductase (4D2) Rabbit Monoclonal Antibody
产品货号: AMRe09715



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 Reductase (4D2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POR
别名 (Alternative Names)	CPR; CYPOR; P450R; por;
基因 ID (Gene ID)	5447.0
蛋白 ID (SwissProt ID)	P16435.A synthetic peptide of human Cytochrome P450 Reductase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:100,ICC/IF 1:200-1:500,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	77kDa

研究背景 (Background)


产品名称: Cytochrome P450 Reductase (4D2) Rabbit Monoclonal Antibody
产品货号: AMRe09715

This enzyme is required for electron transfer from NADP to cytochrome P450 in microsomes. It can also provide electron transfer to heme oxygenase and cytochrome B5. This enzyme is required for electron transfer from NADP to cytochrome P450 in microsomes. It can also provide electron transfer to heme oxygenase and cytochrome B5.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot detection of Cytochrome P450 Reductase in K562, Rat Brain, C6, 3T3 cell lysates using Cytochrome P450 Reductase antibody (1:1000 diluted).

注意事项 (Note)

For research use only.