

产品名称: Cytokeratin 15 (19D17) Rabbit Monoclonal Antibody
产品货号: AMRe09729



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 15 (19D17) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT15
别名 (Alternative Names)	CK-15; CK15; cytokeratin15; K15; K1C15; K1CO; keratin 15; Keratin-15; keratin-15, beta; KRT15; KRTB; type I cytoskeletal 15;
基因 ID (Gene ID)	3866.0
蛋白 ID (SwissProt ID)	P19012.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	49kDa

产品名称: Cytokeratin 15 (19D17) Rabbit Monoclonal Antibody
产品货号: AMRe09729




研究背景 (Background)

The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. There are two types of cytoskeletal and microfibrillar keratin: type I (acidic; 40-55 kDa) [K9 to K20] and type II (neutral to basic; 56-70 kDa) [K1 to K8]. Both a basic and an acidic keratin are required for filament assembly. Note: This description may include information from UniProtKB.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cytokeratin 15 expression in human fetal thymus lysate.

注意事项 (Note)

For research use only .