

产品名称: Cytokeratin 19 (8L9) Rabbit Monoclonal Antibody
产品货号: AMRe09738



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 19 (8L9) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IF-P
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT19
别名 (Alternative Names)	KRT19; K1CS; CK19; K19; CK 19; Cytokeratin 19; MGC15366;
基因 ID (Gene ID)	3880.0
蛋白 ID (SwissProt ID)	P08727 .A synthetic peptide of human Cytokeratin 19

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:5000,IHC 1:20-1:100,ICC/IF 1:20-1:50,IF-P 1:20-1:50
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)

Cytokeratin 19, also known as KRT19, CK19, CK19, K1CS, MGC15366. Entrez Protein NP_002267. It is a member of the


产品名称: Cytokeratin 19 (8L9) Rabbit Monoclonal Antibody
产品货号: AMRe09738



keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Cytokeratin 19 in K562 cell lysates using Cytokeratin 19 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .