

产品名称: Cytokeratin 4 (11K16) Rabbit Monoclonal Antibody
产品货号: AMRe09744



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 4 (11K16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT4
别名 (Alternative Names)	CK 4; CK-4; CYK4; Cytokeratin4; FLJ31692; K4; Keratin 4; Keratin; Keratin type II cytoskeletal 4; Keratin4; Krt4; type II cytoskeletal 4; Type-II keratin Kb4;
基因 ID (Gene ID)	3851.0
蛋白 ID (SwissProt ID)	P19013.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:20-1:50
蛋白分子量 (Molecular Weight)	56kDa

产品名称: Cytokeratin 4 (11K16) Rabbit Monoclonal Antibody
产品货号: AMRe09744




研究背景 (Background)

Cytokeratins are a subfamily of intermediate filament keratins that are characterized by a remarkable biochemical diversity, which is represented in human epithelial tissues by at least 20 different polypeptides.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Cytokeratin 4 expression in A431 cell lysate.

注意事项 (Note)

For research use only .