

产品名称: Cytokeratin 5 (18R6) Rabbit Monoclonal Antibody
产品货号: AMRe09746



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 5 (18R6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT5
别名 (Alternative Names)	58 kDa cytokeratin; CK-5; CK5; Cytokeratin-5; DDD1; EBS2; K2C5; K5; keratin 5; Krt5;
基因 ID (Gene ID)	3852.0
蛋白 ID (SwissProt ID)	P13647.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:100-1:200,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	62kDa

产品名称: Cytokeratin 5 (18R6) Rabbit Monoclonal Antibody
产品货号: AMRe09746


研究背景 (Background)

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in the basal layer of the epidermis with family member KRT14.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Cytokeratin 5 expression in A431 cell lysate.

注意事项 (Note)

For research use only .