

产品名称: DAZAP1 (14U10) Rabbit Monoclonal Antibody  
产品货号: AMRe09799



## 产品概述 (Summary)

|                        |                                           |
|------------------------|-------------------------------------------|
| 产品名称 (Production Name) | DAZAP1 (14U10) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody    |
| 宿主 (Host)              | Rabbit                                    |
| 应用 (Application)       | WB,IP                                     |
| 种属反应性 (Reactivity)     | Human,Rat                                 |

## 产品性能 (Performance)

|                     |                                                                                                                             |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                                                  |
| 形式 (Form)           | Liquid                                                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| 纯化方式 (Purification) | Affinity purification                                                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                            |
|------------------------|--------------------------------------------|
| 基因名 (Gene Name)        | DAZAP1                                     |
| 别名 (Alternative Names) | DAZ associated protein 1; Dazap1;          |
| 基因 ID (Gene ID)        | 26528.0                                    |
| 蛋白 ID (SwissProt ID)   | Q96EP5.Recombinant protein of human DAZAP1 |

## 产品应用 (Application)

|                          |                                |
|--------------------------|--------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:2000-1:10000,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | 43kDa                          |

## 研究背景 (Background)

RNA-binding protein, which may be required during spermatogenesis. RNA-binding protein, which may be required during

产品名称: DAZAP1 (14U10) Rabbit Monoclonal Antibody  
产品货号: AMRe09799




spermatogenesis.

## 研究领域 (Research Area)

Epigenetics and Nuclear Signaling

## 图片 (Image Data)



Western blot detection of DAZAP1 in K562,C6,3T3,HeLa cell lysates using DAZAP1 antibody(1:1000 diluted).

## 注意事项 (Note)

For research use only .