

产品名称: DAZL (11G2) Rabbit Monoclonal Antibody  
产品货号: AMRe09800



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | DAZL (11G2) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,IP                       |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                                                     |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 纯化方式 (Purification) | Affinity purification                                                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                     |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | DAZL                                                                                                                                                                |
| 别名 (Alternative Names) | DAZ homolog; DAZ like autosomal; DAZH; DAZL; DAZL1; DAZLA; Deleted in azoospermia like 1; Germline specific RNA binding protein; SPGY like autosomal; SPGYLA; Tpx2; |
| 基因 ID (Gene ID)        | 1618.0                                                                                                                                                              |
| 蛋白 ID (SwissProt ID)   | Q92904.                                                                                                                                                             |

## 产品应用 (Application)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:1000-1:5000,ICC/IF 1:200-1:500,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | 33kDa                                                             |

产品名称: DAZL (11G2) Rabbit Monoclonal Antibody  
产品货号: AMRe09800




## 研究背景 (Background)

RNA-binding protein, which is essential for gametogenesis. Plays a central role during spermatogenesis. May act by binding to the 3'-UTR of mRNA and thereby regulating the translation of key transcripts. RNA-binding protein, which is essential for gametogenesis in both males and females. Plays a central role during spermatogenesis. Acts by binding to the 3'-UTR of mRNA, specifically recognizing GUU triplets, and thereby regulating the translation of key transcripts (By similarity).

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of DAZL expression in human testis lysate.

## 注意事项 (Note)

For research use only .