

产品名称: Dopamine Receptor D3 (6H18) Rabbit Monoclonal Antibody
产品货号: AMRe10117



产品概述 (Summary)

产品名称 (Production Name)	Dopamine Receptor D3 (6H18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DRD3
别名 (Alternative Names)	D(3) dopamine receptor; Dopamine D3 receptor; DRD3; ETM1; FET1;
基因 ID (Gene ID)	1814.0
蛋白 ID (SwissProt ID)	P35462.A synthetic peptide of human Dopamine Receptor D3/DRD3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)

产品名称: Dopamine Receptor D3 (6H18) Rabbit Monoclonal Antibody
产品货号: AMRe10117




This is one of the five types (D1 to D5) of receptors for dopamine. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase. Promotes cell proliferation. Dopamine receptor whose activity is mediated by G proteins which inhibit adenylyl cyclase. Promotes cell proliferation.

研究领域 (Research Area)

Neuroactive ligand-receptor interaction;

图片 (Image Data)



Western blot analysis of extracts from Mouse brain tissue using RM5657 at 1:1000.

注意事项 (Note)

For research use only .