

产品名称: DPF2 (12O14) Rabbit Monoclonal Antibody
产品货号: AMRe10126



产品概述 (Summary)

产品名称 (Production Name)	DPF2 (12O14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DPF2
别名 (Alternative Names)	BAF45D; DPF2; REQ; UBID4;
基因 ID (Gene ID)	5977.0
蛋白 ID (SwissProt ID)	Q92785.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:100-1:200,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)


产品名称: DPF2 (12O14) Rabbit Monoclonal Antibody
产品货号: AMRe10126



May be a transcription factor required for the apoptosis response following survival factor withdrawal from myeloid cells. Might also have a role in the development and maturation of lymphoid cells. Plays an active role in transcriptional regulation by binding modified histones H3 and H4 (PubMed:28533407, PubMed:27775714). Is a negative regulator of myeloid differentiation of hematopoietic progenitor cells (PubMed:28533407). Might also have a role in the development and maturation of lymphoid cells (By similarity). Involved in the regulation of non-canonical NF-kappa-B pathway (PubMed:20460684).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DPF2 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .