

产品名称: eEF1A1 (19R16) Rabbit Monoclonal Antibody
产品货号: AMRe10308



产品概述 (Summary)

产品名称 (Production Name)	eEF1A1 (19R16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1A1
别名 (Alternative Names)	EEF1A1; CCS-3; CCS3; EEF-1; EEF1A; EF-Tu; EF1A; GRAF-1EF; LENG7; PTI1;eEF1A-1;
基因 ID (Gene ID)	1915.0
蛋白 ID (SwissProt ID)	P68104.A synthetic peptide of human eEF1A1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:200,ICC/IF 1:100-1:200,FC 1:50-1:100,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	50kDa

产品名称: eEF1A1 (19R16) Rabbit Monoclonal Antibody
产品货号: AMRe10308




研究背景 (Background)

This protein promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis. This protein promotes the GTP-dependent binding of aminoacyl- tRNA to the A-site of ribosomes during protein biosynthesis. Plays a role in the positive regulation of IFNG transcription in T-helper 1 cells as part of an IFNG promoter-binding complex with TXK and PARP1 (PubMed: <http://www.uniprot.org/citations/17177976> target="_blank">17177976).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of eEF1A1 in K562,C6,3T3,HeLa cell lysates using eEF1A1 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .