

产品名称: EGFR (ErbB 1) (9C2) Rabbit Monoclonal Antibody  
产品货号: AMRe10340



## 产品概述 (Summary)

|                        |                                                |
|------------------------|------------------------------------------------|
| 产品名称 (Production Name) | EGFR (ErbB 1) (9C2) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody         |
| 宿主 (Host)              | Rabbit                                         |
| 应用 (Application)       | WB,IHC,ICC/IF,FC,IP                            |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                                |

## 产品性能 (Performance)

|                     |                                                                                                                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                                                     |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |
| 纯化方式 (Purification) | Affinity purification                                                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | EGFR                                                                                                     |
| 别名 (Alternative Names) | ERBB1; Epidermal growth factor receptor precursor; Receptor protein-tyrosine kinase ErbB-1; kinase EGFR; |
| 基因 ID (Gene ID)        | 1956.0                                                                                                   |
| 蛋白 ID (SwissProt ID)   | P00533.                                                                                                  |

## 产品应用 (Application)

|                          |                                                                               |
|--------------------------|-------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,IHC 1:50-1:200,ICC/IF 1:200-1:500,FC 1:20-1:200,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | 134kDa                                                                        |

产品名称: EGFR (ErbB 1) (9C2) Rabbit Monoclonal Antibody  
产品货号: AMRe10340

---

## 研究背景 (Background)

EGFR is a receptor tyrosine kinase. Receptor for epidermal growth factor (EGF) and related growth factors including TGF- $\alpha$ , amphiregulin, betacellulin, heparin-binding EGF-like growth factor, GP30 and vaccinia virus growth factor. Is involved in the control of cell growth and differentiation. Receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses (PubMed:<a href="http://www.uniprot.org/citations/2790960" target="\_blank">2790960</a>, PubMed:<a href="http://www.uniprot.org/citations/10805725" target="\_blank">10805725</a>, PubMed:<a href="http://www.uniprot.org/citations/27153536" target="\_blank">27153536</a>). Known ligands include EGF, TGFA/TGF- $\alpha$ , AREG, epigen/EPGN, BTC/betacellulin, epiregulin/EREG and HBEGF/heparin-binding EGF (PubMed:<a href="http://www.uniprot.org/citations/2790960" target="\_blank">2790960</a>, PubMed:<a href="http://www.uniprot.org/citations/7679104" target="\_blank">7679104</a>, PubMed:<a href="http://www.uniprot.org/citations/8144591" target="\_blank">8144591</a>, PubMed:<a href="http://www.uniprot.org/citations/9419975" target="\_blank">9419975</a>, PubMed:<a href="http://www.uniprot.org/citations/15611079" target="\_blank">15611079</a>, PubMed:<a href="http://www.uniprot.org/citations/12297049" target="\_blank">12297049</a>, PubMed:<a href="http://www.uniprot.org/citations/27153536" target="\_blank">27153536</a>, PubMed:<a href="http://www.uniprot.org/citations/20837704" target="\_blank">20837704</a>, PubMed:<a href="http://www.uniprot.org/citations/17909029" target="\_blank">17909029</a>). Ligand binding triggers receptor homo- and/or heterodimerization and autophosphorylation on key cytoplasmic residues. The phosphorylated receptor recruits adapter proteins like GRB2 which in turn activates complex downstream signaling cascades. Activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PLCgamma-PKC and STATs modules (PubMed:<a href="http://www.uniprot.org/citations/27153536" target="\_blank">27153536</a>). May also activate the NF-kappa-B signaling cascade (PubMed:<a href="http://www.uniprot.org/citations/11116146" target="\_blank">11116146</a>). Also directly phosphorylates other proteins like RGS16, activating its GTPase activity and probably coupling the EGF receptor signaling to the G protein-coupled receptor signaling (PubMed:<a href="http://www.uniprot.org/citations/11602604" target="\_blank">11602604</a>). Also phosphorylates MUC1 and increases its interaction with SRC and CTNNB1/beta-catenin (PubMed:<a href="http://www.uniprot.org/citations/11483589" target="\_blank">11483589</a>). Positively regulates cell migration via interaction with CCDC88A/GIV which retains EGFR at the cell membrane following ligand stimulation, promoting EGFR signaling which triggers cell migration (PubMed:<a href="http://www.uniprot.org/citations/20462955" target="\_blank">20462955</a>). Plays a role in enhancing learning and memory performance (By similarity).


## 研究领域 (Research Area)

Signal Transduction

产品名称: EGFR (ErbB 1) (9C2) Rabbit Monoclonal Antibody  
产品货号: AMRe10340



## 图片 (Image Data)



Western blot analysis of EGFR in HeLa cell lysate.

## 注意事项 (Note)

For research use only .