

产品名称: Ezrin (6V7) Rabbit Monoclonal Antibody
产品货号: AMRe10670



产品概述 (Summary)

产品名称 (Production Name)	Ezrin (6V7) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EZR
别名 (Alternative Names)	CVIL; Cytovillin; Villin2; p81; EZR; VIL2;
基因 ID (Gene ID)	7430.0
蛋白 ID (SwissProt ID)	P15311.A synthetic peptide of human Ezrin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:500,IF-P 1:200-1:500
蛋白分子量 (Molecular Weight)	69kDa

研究背景 (Background)

Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for

产品名称: Ezrin (6V7) Rabbit Monoclonal Antibody
产品货号: AMRe10670




the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis. Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of Ezrin in HeLa,A549 cell lysates using Ezrin antibody(1:1000 diluted).

注意事项 (Note)

For research use only .