

产品名称: FADD (18S3) Rabbit Monoclonal Antibody
产品货号: AMRe10793



产品概述 (Summary)

产品名称 (Production Name)	FADD (18S3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	FADD protein; GIG3; Growth inhibiting gene 3 protein; Mediator of receptor induced toxicity; MORT1;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q61160(mouse)(human:Q13158).A synthetic peptide of mouse FADD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	23kDa

研究背景 (Background)


产品名称: FADD (18S3) Rabbit Monoclonal Antibody
产品货号: AMRe10793



Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of FADD in SHSY-5Y cell lysates using FADD antibody(1:1000 diluted).

注意事项 (Note)

For research use only .