

产品名称: Fas (6T3) Rabbit Monoclonal Antibody
产品货号: AMRe10832



产品概述 (Summary)

产品名称 (Production Name)	Fas (6T3) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAS
别名 (Alternative Names)	FAS; ALPS1A; APO-1; APT1; CD95; FAS1; FASTM; TNFRSF6;
基因 ID (Gene ID)	355.0
蛋白 ID (SwissProt ID)	P25445.A synthetic peptide of human Fas

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	38kDa

研究背景 (Background)

Receptor for TNFSF6/FASLG. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-

产品名称: Fas (6T3) Rabbit Monoclonal Antibody
产品货号: AMRe10832




inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis. Receptor for TNFSF6/FASLG. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis. FAS-mediated apoptosis may have a role in the induction of peripheral tolerance, in the antigen-stimulated suicide of mature T-cells, or both. The secreted isoforms 2 to 6 block apoptosis (in vitro).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of Fas in HeLa, Jurkat cell lysates using Fas antibody (1:500 diluted).

注意事项 (Note)

For research use only.