

产品名称: Ferritin Light Chain (10S14) Rabbit Monoclonal Antibody
产品货号: AMRe10906



产品概述 (Summary)

产品名称 (Production Name)	Ferritin Light Chain (10S14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FTL
别名 (Alternative Names)	Ferritin L chain; Ferritin L subunit; Ferritin light chain; Ferritin light polypeptide; ferritin light polypeptide like 3; FTL; LFTD; NBIA 3;
基因 ID (Gene ID)	2512.0
蛋白 ID (SwissProt ID)	P02792.A synthetic peptide of human Ferritin light chain

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	20kDa

产品名称: Ferritin Light Chain (10S14) Rabbit Monoclonal Antibody
产品货号: AMRe10906


研究背景 (Background)

Stores iron in a soluble, non-toxic, readily available form. Important for iron homeostasis. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Also plays a role in delivery of iron to cells. Mediates iron uptake in capsule cells of the developing kidney. Stores iron in a soluble, non-toxic, readily available form. Important for iron homeostasis. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Also plays a role in delivery of iron to cells. Mediates iron uptake in capsule cells of the developing kidney (By similarity).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of extracts from MCF-7 cells using RM5355 at 1:1000.

注意事项 (Note)

For research use only .