

产品名称: fetal hemoglobin (14V10) Rabbit Monoclonal Antibody
产品货号: AMRe10908



产品概述 (Summary)

产品名称 (Production Name)	fetal hemoglobin (14V10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HBG1; HBG2; HBGA; HBGR; Hemoglobin gamma A; Hemoglobin gamma G;
别名 (Alternative Names)	Hemoglobin subunit gamma 1; Hemoglobin subunit gamma 2; HSGGL1; PRO2979;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P69891/P69892.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
----------------------	------------------

产品名称: fetal hemoglobin (14V10) Rabbit Monoclonal Antibody
产品货号: AMRe10908




蛋白分子量 (Molecular Weight) 12kDa

研究背景 (Background)

The gamma globin genes (HBG1 and HBG2) are normally expressed in the fetal liver, spleen and bone marrow. Two gamma chains together with two alpha chains constitute fetal hemoglobin (HbF) which is normally replaced by adult hemoglobin (HbA) at birth.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of fetal hemoglobin expression in K562 cell lysate.

注意事项 (Note)

For research use only .