

产品名称: FHL2 (4I16) Rabbit Monoclonal Antibody
产品货号: AMRe10960



产品概述 (Summary)

产品名称 (Production Name)	FHL2 (4I16) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FHL2
别名 (Alternative Names)	AAG11; DRAL; Fhl2; SLIM3;
基因 ID (Gene ID)	2274.0
蛋白 ID (SwissProt ID)	Q14192.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:500, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	32kDa

研究背景 (Background)


产品名称: FHL2 (4I16) Rabbit Monoclonal Antibody
产品货号: AMRe10960



May function as a molecular transmitter linking various signaling pathways to transcriptional regulation. Negatively regulates the transcriptional repressor E4F1 and may function in cell growth. May function as a molecular transmitter linking various signaling pathways to transcriptional regulation. Negatively regulates the transcriptional repressor E4F1 and may function in cell growth. Inhibits the transcriptional activity of FOXO1 and its apoptotic function by enhancing the interaction of FOXO1 with SIRT1 and FOXO1 deacetylation. Negatively regulates the calcineurin/NFAT signaling pathway in cardiomyocytes (PubMed:28717008).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FHL2 expression in SW480 cell lysate.

注意事项 (Note)

For research use only .