

产品名称: Fos B (18J10) Rabbit Monoclonal Antibody
产品货号: AMRe11069



产品概述 (Summary)

产品名称 (Production Name)	Fos B (18J10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOSB
别名 (Alternative Names)	G0/G1 switch regulatory protein 3, G0S3, Protein fosB; GOSB; Oncogene FOSB; Activator protein 1;
基因 ID (Gene ID)	2354.0
蛋白 ID (SwissProt ID)	P53539.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:200-1:500,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	36kDa

产品名称: Fos B (18J10) Rabbit Monoclonal Antibody
产品货号: AMRe11069




研究背景 (Background)

FosB interacts with Jun proteins enhancing their DNA binding activity. The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. FosB interacts with Jun proteins enhancing their DNA binding activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FosB expression in C6 cell lysate treated with serum.

注意事项 (Note)

For research use only .