

产品名称: Galectin 1 (5D19) Rabbit Monoclonal Antibody
产品货号: AMRe11265



产品概述 (Summary)

产品名称 (Production Name)	Galectin 1 (5D19) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LGALS1
别名 (Alternative Names)	GAL1; Galaptin; Galbp; Galectin; Galectin1; GBP; HBL; HLBP14; HPL; L14; Lect14; LGALS1;
基因 ID (Gene ID)	3956.0
蛋白 ID (SwissProt ID)	P09382.A synthetic peptide of human Galectin 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,FC 1:20-1:50,IP 1:20-1:50,IF-P 1:50-1:100
蛋白分子量 (Molecular Weight)	15kDa

研究背景 (Background)

产品名称: Galectin 1 (5D19) Rabbit Monoclonal Antibody
产品货号: AMRe11265




May regulate apoptosis, cell proliferation and cell differentiation. Binds beta-galactoside and a wide array of complex carbohydrates. Inhibits CD45 protein phosphatase activity and therefore the dephosphorylation of Lyn kinase. Lectin that binds beta-galactoside and a wide array of complex carbohydrates. Plays a role in regulating apoptosis, cell proliferation and cell differentiation. Inhibits CD45 protein phosphatase activity and therefore the dephosphorylation of Lyn kinase. Strong inducer of T-cell apoptosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot detection of Galectin 1 in CHO-K1, C6 cell lysates using Galectin 1 antibody (1:1000 diluted).

注意事项 (Note)

For research use only.