

产品名称: Galectin 3 (3E18) Rabbit Monoclonal Antibody
产品货号: AMRe11268



产品概述 (Summary)

产品名称 (Production Name)	Galectin 3 (3E18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LGALS3
别名 (Alternative Names)	Carbohydrate binding protein 35; CBP35; GAL3; Galactose-specific lectin 3; GALBP; Galectin3; GALIG; GBP; IgE binding protein; L31; LGALS2; LGALS3;
基因 ID (Gene ID)	3958.0
蛋白 ID (SwissProt ID)	P17931.A synthetic peptide of human Galectin 3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:500-1:1000,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	26kDa

研究背景 (Background)


产品名称: Galectin 3 (3E18) Rabbit Monoclonal Antibody
产品货号: AMRe11268

Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration. Together with DMBT1, required for terminal differentiation of columnar epithelial cells during early embryogenesis. Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration. Together with DMBT1, required for terminal differentiation of columnar epithelial cells during early embryogenesis (By similarity). In the nucleus: acts as a pre-mRNA splicing factor. Involved in acute inflammatory responses including neutrophil activation and adhesion, chemoattraction of monocytes macrophages, opsonization of apoptotic neutrophils, and activation of mast cells. Together with TRIM16, coordinates the recognition of membrane damage with mobilization of the core autophagy regulators ATG16L1 and BECN1 in response to damaged endomembranes.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot detection of Galectin 3 in 3T3, HeLa cell lysates using Galectin 3 antibody (1:1000 diluted).

注意事项 (Note)

For research use only .