

产品名称: G-CSF (11N11) Rabbit Monoclonal Antibody
产品货号: AMRe11370



产品概述 (Summary)

产品名称 (Production Name)	G-CSF (11N11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CSF3
别名 (Alternative Names)	CSF beta; CSF3; CSF3OS; Csfg; Filgrastim; G-CSF; GCSA; GCSF; MGI 2; Pluripoietin;
基因 ID (Gene ID)	1440.0
蛋白 ID (SwissProt ID)	P09919.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:200-1:500, FC 1:100-1:200, IP 1:50-1:200
蛋白分子量 (Molecular Weight)	22kDa

产品名称: G-CSF (11N11) Rabbit Monoclonal Antibody
产品货号: AMRe11370




研究背景 (Background)

Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. This CSF induces granulocytes. Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. This CSF induces granulocytes.

研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Western blot analysis of G-CSF expression in K562 cell lysate.

注意事项 (Note)

For research use only .